

Instant Hand Sanitizer Unscented

SDS Preparation Date (mm/dd/yyyy): 10/14/2019

Page 1 of 9

SAFETY DATA SHEET**SECTION 1. IDENTIFICATION****Product identifier used on the label**: **Instant Hand Sanitizer Unscented****Other means of identification** : 9900078
8 fl.oz bottle**Recommended use of the chemical and restrictions on use**: Personal Care Applications
Restriction on use: None known**Chemical family** : Mixture.**Name, address, and telephone number
of the supplier:****Safetec of America, Inc.**887 Kensington Avenue
Buffalo, NY, USA
14215

Supplier's Telephone # : (716) 895 1822

24 Hr. Emergency Tel # : 1-800-255-3924**Name, address, and telephone number of
the manufacturer:**

Refer to supplier

SECTION 2. HAZARDS IDENTIFICATION**Classification of the chemical**

Clear gel liquid. Alcohol odor

This product is sold as a personal care item that is safe for consumers and other users under normal and reasonably foreseeable use. As such, it is not regulated under Hazcom 2012/WHMIS 2015 labeling and SDS requirements do not apply.

Label elements*Hazard pictogram(s)*

None required under U.S. OSHA Hazcom 2012 and Canadian WHMIS 2015 regulations.

Signal Word

None required under OSHA HazCom2012 and 2015 Canadian WHMIS regulations.

Hazard statement(s)

None required under OSHA HazCom2012 and 2015 Canadian WHMIS regulations.

Precautionary statement(s)

None required under OSHA HazCom2012 and 2015 Canadian WHMIS regulations.

Other hazards

None known.

SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS

Mixture

<u>Chemical name</u>	<u>Common name and synonyms</u>	<u>CAS #</u>	<u>Concentration (% by weight)</u>
Ethyl alcohol	Ethyl alcohol Ethyl hydrate	64-17-5	45.0 - 70.0

Instant Hand Sanitizer Unscented

SDS Preparation Date (mm/dd/yyyy): 10/14/2019

Page 2 of 9

SAFETY DATA SHEET

Triethanolamine	N/Av	102-71-6	0.1 - 1.0
-----------------	------	----------	-----------

The exact concentrations of the above listed chemicals are being withheld as a trade secret.

SECTION 4. FIRST-AID MEASURES

Description of first aid measures

- Ingestion* : Harmful effects are not expected under normal usage. In the case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).
- Inhalation* : Harmful effects are not expected under normal usage. In the case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).
- Skin contact* : Harmful effects are not expected under normal usage. If skin irritation or rash occurs: Get medical advice/attention.
- Eye contact* : If in eyes: Rinse cautiously with water for several minutes. If eye irritation persists: get medical advice/attention.

Most important symptoms and effects, both acute and delayed

- : None expected, when used as intended.

Indication of any immediate medical attention and special treatment needed

- : Treat symptomatically.

SECTION 5. FIRE-FIGHTING MEASURES

Extinguishing media

- Suitable extinguishing media* : Dry chemicals. Foam. Carbon dioxide (CO₂). Water or fog spray.
- Unsuitable extinguishing media* : Do not use water jet as extinguisher, as this will spread the fire.

Special hazards arising from the substance or mixture / Conditions of flammability

- : Highly flammable liquid and vapor.

Flammability classification (OSHA 29 CFR 1910.106)

- : Flammable liquids - Category 2

Hazardous combustion products

- : Carbon oxides

Special protective equipment and precautions for firefighters

- Protective equipment for fire-fighters* : Firefighters should wear an approved full-face, self-contained breathing apparatus (SCBA) and impervious clothing.
- Special fire-fighting procedures* : Cool containers exposed to heat with water spray and remove container, if no risk is involved.

SECTION 6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

- : Eliminate all ignition sources (no smoking, flares, sparks, or flames in immediate area). Wear appropriate personal protective equipment.

Environmental precautions

- : Avoid release to the environment.

Methods and material for containment and cleaning up

- : Wipe up with absorbent material (e.g. cloth, fleece).

Special spill response procedures

- : Not applicable.

SECTION 7. HANDLING AND STORAGE

Precautions for safe handling

Instant Hand Sanitizer Unscented

SDS Preparation Date (mm/dd/yyyy): 10/14/2019

Page 3 of 9

SAFETY DATA SHEET

: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources - No smoking. Do not ingest or swallow. Observe good industrial hygiene practices.
Conditions for safe storage : Keep containers closed when not in use.
Incompatible materials : Strong oxidizing agents. Acids.

SECTION 8. EXPOSURE CONTROLS / PERSONAL PROTECTION

<u>Exposure Limits:</u>				
<u>Chemical Name</u>	<u>ACGIH TLV</u>		<u>OSHA PEL</u>	
	<u>TWA</u>	<u>STEL</u>	<u>PEL</u>	<u>STEL</u>
Ethyl alcohol	N/Av	1000 ppm	1000 ppm (1900 mg/m ³)	N/Av
Triethanolamine	5 mg/m ³	N/Av	N/Av	N/Av

Exposure controls

Ventilation and engineering measures

: None required under normal conditions.

Respiratory protection : None required under normal conditions.

Skin protection : None required under normal conditions.

Eye / face protection : None required under normal conditions.

Other protective equipment : None required when used as intended.

General hygiene considerations

: Always observe good personal hygiene measures, such as washing after handling the material and before eating, drinking, and/or smoking.

SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance : Clear colourless liquid.

Odour : Alcohol odour.

Odour threshold : Not available.

pH : 7

Melting Point/Freezing point : Not available.

Initial boiling point and boiling range

: 82.22-100°C (180-212°F)

Flash point : 11.7°C (53°F)

Flashpoint (Method) : Closed cup

Evaporation rate (BuAe = 1) : 2

Flammability (solid, gas) : Not applicable.

Lower flammable limit (% by vol.)

: 2%

Upper flammable limit (% by vol.)

: 12%

Oxidizing properties : Not applicable

Explosive properties : Not explosive

Vapour pressure : Not available.

Vapour density : 2.5

Relative density / Specific gravity

: 0.79

Instant Hand Sanitizer Unscented

SDS Preparation Date (mm/dd/yyyy): 10/14/2019

Page 4 of 9

SAFETY DATA SHEET

Solubility in water : Complete
Other solubility(ies) : Not available.
Partition coefficient: n-octanol/water or Coefficient of water/oil distribution
: Not available.
Auto-ignition temperature : Not available.
Decomposition temperature : Not available.
Viscosity : Not available.
Volatiles (% by weight) : Not available.
Volatile organic Compounds (VOC's)
: Not available.
Absolute pressure of container
: Not available.
Flame projection length : Not available.
Other physical/chemical comments
: None known or reported by the manufacturer.

SECTION 10. STABILITY AND REACTIVITY

Reactivity : The product is stable and non-reactive under normal conditions of use, storage and transport.
Chemical stability : Material is stable under normal conditions.
Possibility of hazardous reactions
: Hazardous polymerization does not occur.
Conditions to avoid : Keep away from heat, sparks and open flames.
Incompatible materials : See Section 7 (Handling and Storage) for further details.
Hazardous decomposition products
: See Section 5 (Fire Fighting Measures).

SECTION 11. TOXICOLOGICAL INFORMATION**Information on likely routes of exposure:**

Routes of entry inhalation : YES
Routes of entry skin & eye : YES
Routes of entry Ingestion : YES
Routes of exposure skin absorption
: YES

Potential Health Effects:**Signs and symptoms of short-term (acute) exposure***Sign and symptoms Inhalation*

: None expected, when used as intended.

Sign and symptoms ingestion

: None expected, when used as intended.

Sign and symptoms skin

: None expected, when used as intended.

Sign and symptoms eyes

: None expected, when used as intended. Direct eye contact may cause irritation. Symptoms may include tearing, redness and discomfort.

Potential Chronic Health Effects

: None expected, when used as intended.

Mutagenicity

: Not expected to be mutagenic.

Instant Hand Sanitizer Unscented

SDS Preparation Date (mm/dd/yyyy): 10/14/2019

Page 5 of 9

SAFETY DATA SHEET

Carcinogenicity : This product is not considered to be a carcinogen by IARC, ACGIH, NTP, or OSHA.

Reproductive effects & Teratogenicity

: This product is not expected to cause reproductive or developmental effects.

Sensitization to material : Not expected to be a skin or respiratory sensitizer.

Specific target organ effects : None expected, when used as intended.

Medical conditions aggravated by overexposure

: None reported by the manufacturer.

Synergistic materials : None reported by the manufacturer.

Toxicological data : See below for individual ingredient acute toxicity data.

<u>Chemical name</u>	LC₅₀(4hr) <u>inh, rat</u>	LD₅₀	
		<u>(Oral, rat)</u>	<u>(Rabbit, dermal)</u>
Ethyl alcohol	> 32 380 ppm (61 mg/L) (vapour)	7060 mg/kg	> 15 800 mg/kg
Triethanolamine	N/Av	6110 mg/kg	> 2000 mg/kg (No mortality)

Other important toxicological hazards

: None known or reported by the manufacturer.

SECTION 12. ECOLOGICAL INFORMATION

Ecotoxicity : The product is not classified as environmentally hazardous. Avoid release to the environment.

Ecotoxicity data:

<u>Ingredients</u>	CAS No	Toxicity to Fish		
		LC₅₀ / 96h	NOEC / 21 day	M Factor
Ethyl alcohol	64-17-5	> 100 mg/L (Fathead minnow)	N/Av	None.
Triethanolamine	102-71-6	11 800 mg/L (Fathead minnow)	N/Av	None.

<u>Ingredients</u>	CAS No	Toxicity to Daphnia		
		EC₅₀ / 48h	NOEC / 21 day	M Factor
Ethyl alcohol	64-17-5	5012 mg/L (Daphnia magna)	N/Av	None.
Triethanolamine	102-71-6	609.88 mg/L [Ceriodaphnia (water flea)]	16 mg/L	None.

<u>Ingredients</u>	CAS No	Toxicity to Algae		
		EC₅₀ / 96h or 72h	NOEC / 96h or 72h	M Factor
Ethyl alcohol	64-17-5	1000 mg/L/96hr (Green algae)	N/Av	None.
Triethanolamine	102-71-6	216 mg/L/72hr (Green algae)	N/Av	None.

Instant Hand Sanitizer Unscented

SDS Preparation Date (mm/dd/yyyy): 10/14/2019

Page 6 of 9

SAFETY DATA SHEET

Persistence and degradability

: Not available.

Bioaccumulation potential

: Not available.

<u>Components</u>	<u>Partition coefficient n-octanol/water (log Kow)</u>	<u>Bioconcentration factor (BCF)</u>
Ethyl alcohol (CAS 64-17-5)	- 0.31	N/Av
Triethanolamine (CAS 102-71-6)	- 1.59	< 3.9 (common carp)

Mobility in soil : Not available.

Other Adverse Environmental effects

: None known or reported by the manufacturer.






SECTION 13. DISPOSAL CONSIDERATIONS

Handling for Disposal : Dispose of in accordance with local regulations.

Methods of Disposal : Dispose of in accordance with local regulations.

RCRA : Not applicable.

SECTION 14. TRANSPORT INFORMATION

Regulatory Information	UN Number	UN proper shipping name	Transport hazard class(es)	Packing Group	Label
TDG	UN1170	ETHANOL	LQ ground	II	
TDG Additional information	May be shipped as LIMITED QUANTITY when transported in quantities no larger than 1 Litre, in packages not exceeding 30 kg gross mass.				
IMDG	UN1170	ETHANOL	LQ ground	II	
IMDG Additional information	May be shipped as LIMITED QUANTITY when transported in quantities no larger than 1 Litre, in packages not exceeding 30 kg gross mass.				
49CFR/DOT	UN1170	ETHANOL	LQ ground	II	
49CFR/DOT Additional information	May be shipped as LIMITED QUANTITY when transported in quantities no larger than 1 Litre, in packages not exceeding 30 kg gross mass.				
ICAO/IATA	UN1170	ETHANOL	3	II	 
ICAO/IATA Additional information	Refer to ICAO/IATA Packing Instruction				

Special precautions for user : Appropriate advice on safety must accompany the package. Keep away from heat and open flames. - No smoking.

Instant Hand Sanitizer Unscented

SDS Preparation Date (mm/dd/yyyy): 10/14/2019

Page 7 of 9

SAFETY DATA SHEET

Environmental hazards : This product does not meet the criteria for an environmentally hazardous mixture, according to the IMDG Code. See Section 12 for more environmental information.

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code
: Not established.

SECTION 15 - REGULATORY INFORMATION

US Federal Information:

Components listed below are present on the following U.S. Federal chemical lists:

<u>Ingredients</u>	CAS #	TSCA Inventory	CERCLA Reportable Quantity(RQ) (40 CFR 117.302):	SARA TITLE III: Sec. 302, Extremely Hazardous Substance, 40 CFR 355:	SARA TITLE III: Sec. 313, 40 CFR 372, Specific Toxic Chemical	
					Toxic Chemical	de minimus Concentration
Ethyl alcohol	64-17-5	Yes	None.	None.	No	N/Ap
Triethanolamine	102-71-6	Yes	N/Ap	N/Ap	No	N/Ap

SARA TITLE III: Sec. 313, Toxic Chemicals Notification, 40 CFR 372: None. Under SARA Sections 311 and 312, the EPA has established threshold quantities for the reporting of hazardous chemicals. The current thresholds are 500 pounds or the threshold planning quantity (TPQ), whichever is lower, for extremely hazardous substances and 10,000 pounds for all other hazardous chemicals.

US State Right to Know Laws:

The following chemicals are specifically listed by individual States:

<u>Ingredients</u>	CAS #	California Proposition 65		State "Right to Know" Lists					
		Listed	Type of Toxicity	CA	MA	MN	NJ	PA	RI
Ethyl alcohol	64-17-5	No	N/Ap	Yes	Yes	Yes	Yes	Yes	Yes
Triethanolamine	102-71-6	No	N/Ap	No	Yes	Yes	Yes	Yes	Yes

Canadian Information:

All components of this product are on the Canadian DSL list.
WHMIS Classification: See Section 2.

International Information:

Components listed below are present on the following International Inventory list:

<u>Ingredients</u>	CAS #	European EINECS	Australia AICS	Philippines PICCS	Japan ENCS	Korea KECI/KECL	China IECSC	NewZealand IOC
Ethyl alcohol	64-17-5	200-578-6	Present	Present	(2)-202	KE-13217	Present	HSR001144
Triethanolamine	102-71-6	203-049-8	Present	Present	(2)-308	KE-25940	Present	HSR002785

SECTION 16. OTHER INFORMATION

Legend : ACGIH: American Conference of Governmental Industrial Hygienists
CA: California
CAS: Chemical Abstract Services

Instant Hand Sanitizer Unscented**SDS Preparation Date (mm/dd/yyyy): 10/14/2019**

Page 8 of 9

SAFETY DATA SHEET

CEPA: Canadian Environmental Protection Act
CERCLA: Comprehensive Environmental Response, Compensation, and Liability Act of 1980
CFR: Code of Federal Regulations
DOT: Department of Transportation
DSL: Domestic Substances List
EPA: Environmental Protection Agency
HPA: Hazardous Products Act
HSDB: Hazardous Substances Data Bank
IARC: International Agency for Research on Cancer
IATA: International Air Transport Association
ICAO: International Civil Aviation Organisation
IMDG: International Maritime Dangerous Goods
Inh: Inhalation
LC: Lethal Concentration
LD: Lethal Dose
MA: Massachusetts
MN: Minnesota
N/Ap: Not Applicable
N/Av: Not Available
NIOSH: National Institute of Occupational Safety and Health
NJ: New Jersey
NOEC: No observable effect concentration
NTP: National Toxicology Program
OECD: Organisation for Economic Co-operation and Development
OEL: National occupational exposure limits
OSHA: Occupational Safety and Health Administration
PA: Pennsylvania
PEL: Permissible exposure limit
PPE: Personal Protective Equipment
RCRA: Resource Conservation and Recovery Act
RI: Rhode Island
RQ: Reportable Quantity
RTECS: Registry of Toxic Effects of Chemical Substances
SARA: Superfund Amendments and Reauthorization Act
SDS: Safety Data Sheet
STEL: Short Term Exposure Limit
TDG: Canadian Transportation of Dangerous Goods Act & Regulations
TLV: Threshold Limit Values
TWA: Time Weighted Average
WHMIS: Workplace Hazardous Materials Identification System

References

- : 1. ACGIH, Threshold Limit Values for Chemical Substances and Physical Agents & Biological Exposure Indices for 2018.
- 2. International Agency for Research on Cancer Monographs, searched 2019.
- 3. Canadian Centre for Occupational Health and Safety, CCInfoWeb databases, 2019 (Chempendium, HSDB and RTECs).
- 4. Safety Data Sheets from manufacturer.
- 5. US EPA Title III List of Lists - June 2019 version.
- 6. California Proposition 65 List - September 2019 version.
- 7. OECD - The Global Portal to Information on Chemical Substances - eChemPortal, 2019.

Preparation Date (mm/dd/yyyy)

: 10/14/2019

Other special considerations for handling

: Provide adequate information, instruction and training for operators.

Instant Hand Sanitizer Unscented**SDS Preparation Date (mm/dd/yyyy): 10/14/2019**

Page 9 of 9

SAFETY DATA SHEET

<u>Prepared for:</u> Safetec of America Inc. 887 Kensington Avenue Buffalo, NY 14215 Telephone: (716) 895 1822 www.safetec.com	
<u>Prepared by:</u> ICC The Compliance Center Inc. Telephone: (888) 442-9628 (U.S.); (888) 977-4834 (Canada) http://www.thecompliancecenter.com	

DISCLAIMER

This Safety Data Sheet was prepared by ICC The Compliance Center Inc. using information provided by Safetec of America Inc. and CCOHS' Web Information Service. The information in the Safety Data Sheet is offered for your consideration and guidance when exposed to this product. ICC The Compliance Center Inc and Safetec of America Inc. expressly disclaim all expressed or implied warranties and assume no responsibilities for the accuracy or completeness of the data contained herein. The data in this Safety Data Sheet does not apply to use with any other product or in any other process. This Safety Data Sheet may not be changed, or altered in any way without the expressed knowledge and permission of ICC The Compliance Center Inc. and Safetec of America Inc.

END OF DOCUMENT